



# NSUN3 Antibody

<b>Product Code</b>	CSB-PA867156LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9H649
<b>Immunogen</b>	Recombinant Human tRNA (cytosine(34)-C(5))-methyltransferase, mitochondrial protein (1-200AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	tRNA (cytosine(34)-C(5))-methyltransferase, mitochondrial (EC 2.1.1.-) (NOL1/NOP2/Sun domain family member 3), NSUN3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	NSUN3