



NSUN3 Antibody, Biotin conjugated

Product Code	CSB-PA867156LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H649
Immunogen	Recombinant Human tRNA (cytosine(34)-C(5))-methyltransferase, mitochondrial protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	tRNA (cytosine(34)-C(5))-methyltransferase, mitochondrial (EC 2.1.1.-) (NOL1/NOP2/Sun domain family member 3), NSUN3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NSUN3