



RHOF Antibody, FITC conjugated

Product Code	CSB-PA867192LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBH0
Immunogen	Recombinant Human Rho-related GTP-binding protein RhoF protein (1-170AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-related GTP-binding protein RhoF (Rho family GTPase Rif) (Rho in filopodia), RHOF, ARHF RIF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RHOF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.