



PLAC1 Antibody, FITC conjugated

Product Code	CSB-PA867193LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HBJ0
Immunogen	Recombinant Human Placenta-specific protein 1 protein (23-212AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	May play a role in placental development.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Placenta-specific protein 1, PLAC1
Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	PLAC1