





TRIM36 Antibody

Product Code	CSB-PA868285LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQ86
Immunogen	Recombinant Human E3 ubiquitin-protein ligase TRIM36 protein (70-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase TRIM36 (EC 2.3.2.27) (RING finger protein 98) (RING-type E3 ubiquitin transferase TRIM36) (Tripartite motif-containing protein 36) (Zinc-binding protein Rbcc728), TRIM36, RBCC728 RNF98
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM36

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA868285LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA868285LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.