



TRIM36 Antibody, Biotin conjugated

Product Code	CSB-PA868285LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQ86
Immunogen	Recombinant Human E3 ubiquitin-protein ligase TRIM36 protein (70-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase TRIM36 (EC 2.3.2.27) (RING finger protein 98) (RING-type E3 ubiquitin transferase TRIM36) (Tripartite motif-containing protein 36) (Zinc-binding protein Rbcc728), TRIM36, RBCC728 RNF98
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM36
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.