




SACM1L Antibody

Product Code	CSB-PA868317LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NTJ5
Immunogen	Recombinant Human Phosphatidylinositide phosphatase SAC1 protein (301-520AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositide phosphatase SAC1 (EC 3.1.3.-) (Suppressor of actin mutations 1-like protein), SACM1L, KIAA0851
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SACM1L

Image



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA868317LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.