





PLEKHJ1 Antibody

Product Code	CSB-PA868334LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NW61
Immunogen	Recombinant Human Pleckstrin homology domain-containing family J member 1 protein (1-149AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pleckstrin homology domain-containing family J member 1 (PH domain-containing family J member 1) (Guanine nucleotide-releasing protein x), PLEKHJ1, GNRPX
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PLEKHJ1

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA868334LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA868334LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.