



# PLEKHJ1 Antibody, FITC conjugated

|                     |   |
|---------------------|---|
| Product Code        | CSB-PA868334LC01HU  |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| Uniprot No.         | Q9NW61  |
| Immunogen           | Recombinant Human Pleckstrin homology domain-containing family J member 1 protein (1-149AA)   |
| Raised In           | Rabbit  |
| Species Reactivity  | Human   |
| Form                | Liquid  |
| Conjugate           | FITC  |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4  |
| Purification Method | >95%, Protein G purified  |
| Isotype             | IgG   |
| Clonality           | Polyclonal  |
| Alias               | Pleckstrin homology domain-containing family J member 1 (PH domain-containing family J member 1) (Guanine nucleotide-releasing protein x), PLEKHJ1, GNRPX |
| Immunogen Species   | Homo sapiens (Human)  |
| Research Area       | Signal Transduction   |
| Target Names        | PLEKHJ1   |