



OXSM Antibody, Biotin conjugated

Product Code	CSB-PA868340LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NWU1
Immunogen	Recombinant Human 3-oxoacyl-[acyl-carrier-protein] synthase, mitochondrial protein (78-345AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-oxoacyl-[acyl-carrier-protein] synthase, mitochondrial (EC 2.3.1.41) (Beta-ketoacyl-ACP synthase), OXSM
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	OXSM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.