




TAS2R14 Antibody

Product Code	CSB-PA868362LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYV8
Immunogen	Recombinant Human Taste receptor type 2 member 14 protein (283-317AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Taste receptor type 2 member 14 (T2R14) (Taste receptor family B member 1) (TRB1), TAS2R14
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TAS2R14

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA868362LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.