



MRPL27 Antibody, HRP conjugated

Product Code	CSB-PA868378LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0M9
Immunogen	Recombinant Human 39S ribosomal protein L27, mitochondrial protein (1-148AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L27, mitochondrial (L27mt) (MRP-L27) (Mitochondrial large ribosomal subunit protein bL27m), MRPL27
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MRPL27
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.