



MBD5 Antibody

| | |
|---------------------|---|
| Product Code | CSB-PA868388LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9P267 |
| Immunogen | Recombinant Human Methyl-CpG-binding domain protein 5 protein (623-840AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Methyl-CpG-binding domain protein 5 (Methyl-CpG-binding protein MBD5), MBD5, KIAA1461 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | MBD5 |