



MBD5 Antibody, HRP conjugated

Product Code	CSB-PA868388LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P267
Immunogen	Recombinant Human Methyl-CpG-binding domain protein 5 protein (623-840AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methyl-CpG-binding domain protein 5 (Methyl-CpG-binding protein MBD5), MBD5, KIAA1461
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MBD5