





XPR1 Antibody

Product Code	CSB-PA871381LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBH6
Immunogen	Recombinant Human Xenotropic and polytropic retrovirus receptor 1 protein (529-696AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Xenotropic and polytropic retrovirus receptor 1 (Protein SYG1 homolog) (Xenotropic and polytropic murine leukemia virus receptor X3) (X-receptor), XPR1, SYG1 XR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	XPR1


Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA871381LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA871381LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA871381LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)