



HESX1 Antibody

Product Code	CSB-PA871390LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBX0
Immunogen	Recombinant Human Homeobox expressed in ES cells 1 protein (1-185AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox expressed in ES cells 1 (Homeobox protein ANF) (hAnf), HESX1, HANF
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HESX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.