



HESX1 Antibody, HRP conjugated

| | |
|----------------------------|----------------------------------------------------------------------------------|
| Product Code | CSB-PA871390LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UBX0 |
| Immunogen | Recombinant Human Homeobox expressed in ES cells 1 protein (1-185AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Homeobox expressed in ES cells 1 (Homeobox protein ANF) (hAnf), HESX1, HANF |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | HESX1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |