



NAA38 Antibody, Biotin conjugated

Product Code	CSB-PA871557LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRA0
Immunogen	Recombinant Human N-alpha-acetyltransferase 38, NatC auxiliary subunit protein (2-173AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-alpha-acetyltransferase 38, NatC auxiliary subunit (LSM domain-containing protein 1) (Phosphonoformate immuno-associated protein 2), NAA38, LSMD1 MAK31 PFAAP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NAA38
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.