



NAA38 Antibody, Biotin conjugated

Product Code	CSB-PA871557LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRA0
Immunogen	Recombinant Human N-alpha-acetyltransferase 38, NatC auxiliary subunit protein (2-173AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-alpha-acetyltransferase 38, NatC auxiliary subunit (LSM domain-containing protein 1) (Phosphonoformate immuno-associated protein 2), NAA38, LSMD1 MAK31 PFAAP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NAA38