



GINS4 Antibody, HRP conjugated

Product Code	CSB-PA871563LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRT9
Immunogen	Recombinant Human DNA replication complex GINS protein SLD5 protein (2-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA replication complex GINS protein SLD5 (GINS complex subunit 4) [Cleaved into: DNA replication complex GINS protein SLD5, N-terminally processed], GINS4, SLD5
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GINs4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.