





RHPN1-AS1 Antibody

Product Code	CSB-PA871590LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWJ2
Immunogen	Recombinant Human Putative uncharacterized protein encoded by RHPN1-AS1 protein (1-59AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative uncharacterized protein encoded by RHPN1-AS1 (RHPN1 antisense RNA 1) (RHPN1 antisense gene protein 1), RHPN1-AS1, C8orf51
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RHPN1-AS1

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA871590LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA871590LA01HU at dilution of 1:100

