



RHPN1-AS1 Antibody, FITC conjugated

Product Code	CSB-PA871590LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWJ2
Immunogen	Recombinant Human Putative uncharacterized protein encoded by RHPN1-AS1 protein (1-59AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative uncharacterized protein encoded by RHPN1-AS1 (RHPN1 antisense RNA 1) (RHPN1 antisense gene protein 1), RHPN1-AS1, C8orf51
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RHPN1-AS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.