



RHPN1-AS1 Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA871590LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9BWJ2 |
| Immunogen | Recombinant Human Putative uncharacterized protein encoded by RHPN1-AS1 protein (1-59AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Putative uncharacterized protein encoded by RHPN1-AS1 (RHPN1 antisense RNA 1) (RHPN1 antisense gene protein 1), RHPN1-AS1, C8orf51 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | RHPN1-AS1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |