




SLC46A2 Antibody

Product Code	CSB-PA871607LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BY10
Immunogen	Recombinant Human Thymic stromal cotransporter homolog protein (45-78AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IF; Recommended dilution: WB:1:1000-1:5000, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thymic stromal cotransporter homolog (Solute carrier family 46 member 2), SLC46A2, TSCOT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC46A2

Image



Western blot

All lanes: SLC46A2 antibody at 8μg/ml

Lane 1: HT29 whole cell lysate

Lane 2: Mouse liver tissue


Lane 3: Mouse kidney tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 52 kDa

Observed band size: 52 kDa



Immunofluorescent analysis of HepG2 cells using CSB-PA871607LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.