




KRTAP3-2 Antibody

Product Code	CSB-PA871617LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYR7
Immunogen	Recombinant Human Keratin-associated protein 3-2 protein (2-98AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin-associated protein 3-2 (High sulfur keratin-associated protein 3.2) (Keratin-associated protein 3.2), KRTAP3-2, KAP3.2 KRTAP3.2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRTAP3-2

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA871617LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.