



KRTAP3-2 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA871617LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9BYR7 |
| Immunogen | Recombinant Human Keratin-associated protein 3-2 protein (2-98AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Keratin-associated protein 3-2 (High sulfur keratin-associated protein 3.2) (Keratin-associated protein 3.2), KRTAP3-2, KAP3.2 KRTAP3.2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | KRTAP3-2 |