



C22orf23 Antibody

Product Code	CSB-PA871624LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZE7
Immunogen	Recombinant Human UPF0193 protein EVG1 protein (1-217AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UPF0193 protein EVG1, C22orf23
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C22orf23
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.