



NLRP1 Antibody, HRP conjugated

Product Code	CSB-PA871630EB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AM0
Immunogen	Recombinant Human NACHT, LRR and PYD domains-containing protein 1 protein (1-146AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NLRP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.