



NLRP1 Antibody, FITC conjugated

Product Code	CSB-PA871630EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C000
Immunogen	Recombinant Human NACHT, LRR and PYD domains-containing protein 1 protein (1-146AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NACHT, LRR and PYD domains-containing protein 1 (Caspase recruitment domain-containing protein 7) (Death effector filament-forming ced-4-like apoptosis protein) (Nucleotide-binding domain and caspase recruitment domain), NLRP1, CARD7 DEFCAP KIAA0926 NAC NALP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NLRP1