



NBPF3 Antibody, Biotin conjugated

Product Code	CSB-PA872422LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H094
Immunogen	Recombinant Human Neuroblastoma breakpoint family member 3 protein (1-252AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuroblastoma breakpoint family member 3 (Protein AE2) (Protein SHIIIa4), NBPF3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NBPF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.