



R3HDML Antibody

Product Code	CSB-PA872491LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3Y0
Immunogen	Recombinant Human Peptidase inhibitor R3HDML protein (90-253AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidase inhibitor R3HDML (Cysteine-rich secretory protein R3HDML), R3HDML
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	R3HDML
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.