




RHOJ Antibody

Product Code	CSB-PA872495LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4E5
Immunogen	Recombinant Human Rho-related GTP-binding protein RhoJ protein (1-211AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-related GTP-binding protein RhoJ (Ras-like protein family member 7B) (Tc10-like GTP-binding protein), RHOJ, ARHJ RASL7B RHOI TCL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RHOJ

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA872495LA01HU at dilution of 1:100