



RHOJ Antibody, FITC conjugated

Product Code	CSB-PA872495LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4E5
Immunogen	Recombinant Human Rho-related GTP-binding protein RhoJ protein (1-211AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-related GTP-binding protein RhoJ (Ras-like protein family member 7B) (Tc10-like GTP-binding protein), RHOJ, ARHJ RASL7B RHOI TCL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RHOJ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.