



# FBXO44 Antibody, Biotin conjugated

Product Code	CSB-PA872497LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4M3
Immunogen	Recombinant Human F-box only protein 44 protein (1-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box only protein 44 (F-box protein FBX30) (F-box/G-domain protein 3), FBXO44, FBG3 FBX30 FBX44 FBX6A FBXO6A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXO44