



FBXO44 Antibody, Biotin conjugated

Product Code	CSB-PA872497LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4M3
Immunogen	Recombinant Human F-box only protein 44 protein (1-224AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box only protein 44 (F-box protein FBX30) (F-box/G-domain protein 3), FBXO44, FBG3 FBX30 FBX44 FBX6A FBXO6A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXO44
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.