



FNDC4 Antibody, Biotin conjugated

Product Code	CSB-PA872504LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H6D8
Immunogen	Recombinant Human Fibronectin type III domain-containing protein 4 protein (45-167AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibronectin type III domain-containing protein 4 (Fibronectin type III repeat-containing protein 1), FNDC4, FRCP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FNDC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.