



LRRC19 Antibody, Biotin conjugated

Product Code	CSB-PA872509LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H756
Immunogen	Recombinant Human Leucine-rich repeat-containing protein 19 protein (25-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat-containing protein 19, LRRC19
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LRRC19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.