




UBAP1 Antibody

Product Code	CSB-PA873701LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZ09
Immunogen	Recombinant Human Ubiquitin-associated protein 1 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-associated protein 1 (UBAP-1) (Nasopharyngeal carcinoma-associated gene 20 protein), UBAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	UBAP1

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA873701LA01HU at dilution of 1:100