



UBAP1 Antibody, HRP conjugated

Product Code	CSB-PA873701LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZ09
Immunogen	Recombinant Human Ubiquitin-associated protein 1 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-associated protein 1 (UBAP-1) (Nasopharyngeal carcinoma-associated gene 20 protein), UBAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	UBAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.