



NUP85 Antibody, HRP conjugated

Product Code	CSB-PA874821LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BW27
Immunogen	Recombinant Human Nuclear pore complex protein Nup85 protein (201-500AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear pore complex protein Nup85 (85 kDa nucleoporin) (FROUNT) (Nucleoporin Nup75) (Nucleoporin Nup85) (Pericentrin-1), NUP85, NUP75 PCNT1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	NUP85
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.