



NUP85 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA874821LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9BW27 |
| Immunogen | Recombinant Human Nuclear pore complex protein Nup85 protein (201-500AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Nuclear pore complex protein Nup85 (85 kDa nucleoporin) (FROUNT) (Nucleoporin Nup75) (Nucleoporin Nup85) (Pericentrin-1), NUP85, NUP75 PCNT1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Tags & Cell Markers |
| Target Names | NUP85 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |