




C1QTNF6 Antibody

Product Code	CSB-PA874831EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXI9
Immunogen	Recombinant Human Complement C1q tumor necrosis factor-related protein 6 protein (47-278AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement C1q tumor necrosis factor-related protein 6, C1QTNF6, CTRP6
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	C1QTNF6

Image



Western blot


All lanes: C1QTNF6 antibody at 8µg/ml + Mouse brain tissue

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 29, 31, 18 kDa

Observed band size: 29 kDa



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA874831EA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA874831EA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.