



C1QTNF6 Antibody, Biotin conjugated

| | |
|---------------------|---|
| Product Code | CSB-PA874831ED01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9BXI9 |
| Immunogen | Recombinant Human Complement C1q tumor necrosis factor-related protein 6 protein (47-278AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Complement C1q tumor necrosis factor-related protein 6, C1QTNF6, CTRP6 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | C1QTNF6 |