



KRTAP3-3 Antibody

Product Code	CSB-PA874846LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYR6
Immunogen	Recombinant Human Keratin-associated protein 3-3 protein (1-98AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin-associated protein 3-3 (High sulfur keratin-associated protein 3.3) (Keratin-associated protein 3.3), KRTAP3-3, KAP3.3 KRTAP3.3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KRTAP3-3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.