



KRTAP3-3 Antibody, Biotin conjugated

| | |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA874846LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9BYR6 |
| Immunogen | Recombinant Human Keratin-associated protein 3-3 protein (1-98AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Keratin-associated protein 3-3 (High sulfur keratin-associated protein 3.3) (Keratin-associated protein 3.3), KRTAP3-3, KAP3.3 KRTAP3.3 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | KRTAP3-3 |