



OR2H1 Antibody, Biotin conjugated

Product Code	CSB-PA875629LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZK4
Immunogen	Recombinant Human Olfactory receptor 2H1 protein (157-193AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Olfactory receptor 2H1 (Hs6M1-16) (OLFR42A-9004.14/9026.2) (Olfactory receptor 2H6) (Olfactory receptor 2H8) (Olfactory receptor 6-2) (OR6-2) (Olfactory receptor OR6-32), OR2H1, OR2H6 OR2H8
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	OR2H1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.