



TINAGL1 Antibody, Biotin conjugated

Product Code	CSB-PA875630LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZM7
Immunogen	Recombinant Human Tubulointerstitial nephritis antigen-like protein (250-467AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulointerstitial nephritis antigen-like (Glucocorticoid-inducible protein 5) (Oxidized LDL-responsive gene 2 protein) (OLRG-2) (Tubulointerstitial nephritis antigen-related protein) (TIN Ag-related protein) (TIN-Ag-RP), TINAGL1, GIS5 LCN7 OLRG2 TINAGL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TINAGL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.