



SLC6A15 Antibody, HRP conjugated

Product Code	CSB-PA875657LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2J7
Immunogen	Recombinant Human Sodium-dependent neutral amino acid transporter B(0)AT2 protein (357-452AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-dependent neutral amino acid transporter B(0)AT2 (Sodium- and chloride-dependent neurotransmitter transporter NTT73) (Sodium-coupled branched-chain amino-acid transporter 1) (Solute carrier family 6 member 15) (Transporter v7-3), SLC6A15, B0AT2 NTT73 SBAT1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC6A15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.