



SH3GLB1 Antibody, Biotin conjugated

Product Code	CSB-PA878084LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y371
Immunogen	Recombinant Human Endophilin-B1 protein (1-365AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endophilin-B1 (Bax-interacting factor 1) (Bif-1) (SH3 domain-containing GRB2-like protein B1), SH3GLB1, KIAA0491
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SH3GLB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.