




UBE2T Antibody

Product Code	CSB-PA878841LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPD8
Immunogen	Recombinant Human Ubiquitin-conjugating enzyme E2 T protein (1-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-conjugating enzyme E2 T (EC 2.3.2.23) (Cell proliferation-inducing gene 50 protein) (E2 ubiquitin-conjugating enzyme T) (Ubiquitin carrier protein T) (Ubiquitin-protein ligase T), UBE2T
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	UBE2T

Image



Western Blot

All lanes: UBE2T antibody at 10µg/ml

Lane 1: HepG2 whole cell lysate


Lane 2: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 23 kDa

Observed band size: 23 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA878841LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.